


CPAP MASK



Features

1. Can be apply in Emergency, Anaesthesia, Intensive Care, Sleep Apnoea and Renal Failure
2. The soft anatomically shaped cushion will comfortably fit a wide range of patient faces
3. Twin-port design to ensure the patient inhales through the upper and exhales through the lower port
4. Built-in monitoring port to allow connection of supplemental oxygen or pressure monitoring



Size	Description	Top Shell Material	Cushion Material	Valve	Connector
3	Small adult	PC	PVC	Silicone	Inhale 19mm ID/ 22 mm OD
5	Large adult	PC	PVC	Silicone	Exhale 30 mm OD
6	Extra Large adult	PC	PVC	Silicone	

Shelf life 5 years

Option Harness
Size: Adult (L)/Child(M)
Material: silicone (transparent green, transparent white), Rubber (Black)



PEEP valve:
0-10cmH2O/5-20cmH2O
Material: PC housing, silicone valve



Package PE bag per piece

Regulatory ISO5356-1, BS EN 13544-2, ISO-17510-1 ISO-10993-1, ISO-10993-5, ISO-10993-10, EN980